

ASTIGMATIC

SUCCESS

from chair

TO WEAR

Achieve first fit success with  
Bausch + Lomb ULTRA® for Astigmatism



**BAUSCH + LOMB**  
See better. Live better.

BAUSCH + LOMB ULTRA® FOR ASTIGMATISM

# FIT WITH CONFIDENCE


PATIENT INSIGHT

**96%** of patients say their current lenses **don't address low light and glare**<sup>1</sup>



**Bausch + Lomb ULTRA® delivers consistently clear vision**<sup>2</sup>

Designed with **spherical aberration control** optimized at every 0.25D step in both axes to help reduce halos and glare<sup>2,3</sup>



of patients agree this lens **helps reduce halos and glare, even in low light**<sup>4</sup>

# OpticAlign™ Design FOR FIRST FIT SUCCESS

Rotation at  
dispensing of  $\leq 5^\circ$   
for **94%** of patients<sup>2</sup>

## SUCCESS IN THE CHAIR

ECPs rated speed of fitting as **“excellent”**  
or **“very good”** for **92%** of patients<sup>2</sup>

**Provides ballasting in the lower hemisphere  
for a stable fit**

Designed to stabilize with the eye's natural blink movement

## EXPLORE THE HUMAN BLINK



Point your  
camera at the QR  
code to learn more.\*

\*Data plan and roaming charges may apply.

## SUCCESS DURING WEAR

Patients wearing Bausch + Lomb ULTRA® for Astigmatism agree<sup>4</sup>:

91%

feels **comfortable**  
throughout the day

96%

provides **comfortable vision**,  
even when spending hours  
on digital devices

95%

delivers **clear vision**  
when driving at night

# SUCCESSFULLY FIT

even more of your

# ASTIGMATIC PATIENTS

The **only** monthly toric lens with a **-2.75D Cylinder** available as a standard offering

Not made to order | No price differential | Same-day fitting

Around the clock in 10° steps in **all 5 Cylinders**

AXIS	LOW MINUS 0.00D TO -6.00D (IN 0.25D STEPS)					HIGH MINUS -6.50D TO -9.00D (IN 0.50D STEPS)					LOW PLUS +0.25D TO +6.00D (IN 0.25D STEPS)				
	-0.75	-1.25	-1.75	-2.25	-2.75	-0.75	-1.25	-1.75	-2.25	-2.75	-0.75	-1.25	-1.75	-2.25	-2.75
10°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
20°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
30°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
40°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
50°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
60°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
70°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
80°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
90°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
100°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
110°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
120°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
130°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
140°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
150°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
160°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
170°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
180°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

**REFERENCES:** 1. The Multi-sponsor Survey's 2015 Gallup study of the US consumer contact lens market. March 2015. 2. Results from a 7-investigator, multi-site 2-week study of Bausch + Lomb ULTRA® for Astigmatism contact lenses on 157 current soft contact lens wearers. 3. Data on file. Bausch & Lomb Incorporated. Rochester, NY. 4. Results of an online survey with patients who completed an evaluation program for Bausch + Lomb ULTRA® for Astigmatism contact lenses and wore their trial lenses for ≥4 days (n=426). Survey results include patients who strongly agreed, agreed, or slightly agreed (on a 6-point agreement scale) with the surveyed statement, with a margin of error ±2.7%.

®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates.  
©2018 Bausch & Lomb Incorporated. 160-8946E TP15292b